



SL:

- 1) 1 St. Gehäuse Cito-QC, kpl. (550...)
- 2) 1 St. NFP Pr. , Geräteeinbau (205NFP020-G)
- 3) 2 St. Flachdichtung NFP,Standard (501001)
- 4) 2 St. Doppelnippel cito 3/8x3/8, MS (510030)
- 5) 2 St. O-Ring cito in/out, 15x2,5 (502051)
- 6) 1 St. BG Eing.schlauch Vario (577031)
- 7) 1 St. BG Ausg.schlauch Vario (577032)
- 8) 1 St. Vario WS8, Geräteeinbau (120VAR-WS8-G)
- 15) 1 St. Bedienanleitung cito (545...)
- 16) 1 St. VK VHP/Uno inox/cito.i.Griff (610076)
- 17) 1 St. Kart.eins. cito für 610076 (610...)
- 18) 1 St. Kart.überstülper cito (605...)
- 19) 1 St. BG Beipack NFP...-9 (545051)
- 20) 1 St. Patr.kart.NFP/IFP neutral, kpl. (610040)

1	erst.		2.11.15	Hanf

cito-QC Classic

122QC-BS1

